

medartis®

PRECISION IN FIXATION

PRODUCT INFORMATION

MODUS® 2
90°
Screwdriver

MODUS®





Contents


4	Osteosynthesis from a Different Angle
4	Areas of Use
5	The 90° Screwdriver
5	Features
7	Storage
9	Ordering Information

For further information on the MODUS product line, please visit www.medartis.com.

Osteosynthesis from a Different Angle

Areas of Use


- Osteosynthesis at the mandible angle, the ascending ramus and condyle
- Orthognathics (sagittal split/ramus osteotomy)
- Intraoral plate and screw removal




The 90° Screwdriver

Features


- Head height of 8 mm and slim shape of the shaft provide the best possible freedom of movement and overview in situ.
- Standardized INTRA coupling (ISO 3964) allows connection to power systems from different manufacturers, enabling powered drilling.
- For easier manual insertion of the screws, the 90° screwdriver has a 1.75 : 1 transmission.
- Removable holding fork for plates and screws allows one-handed placement and positioning of plates. After use, the holder can be retracted and turned to the side to ensure the best possible visibility in situ.
- Easy to disassemble for automated cleaning and sterilization







Storage

- User-friendly storage of instruments
- System components
 - 1 screwdriver complete with rotation knob and holding fork for screws and plates
 - 1 screwdriver for use as drill
 - Removable storage module for blades and twist drills with color coding





Removable storage module for blades and drills

Ordering Information

90° Screwdriver



M-2440



M-2441

Art. No.	Description	Length	Pieces/Pkg
M-2440	screwdriver complete with screwholding fork M-2442 and rotation knob M-2443	320 mm	1
M-2441		253 mm	1

Handles



M-2443



M-2438

Art. No.	Description	Length	Pieces/Pkg
M-2443	rotation knob	67 mm	1
M-2438	torque limiting knob	58 mm	1

Screwholding Fork



Art. No.	Description	Length	Pieces/Pkg
M-2442	removable	128 mm	1

Screwdriver Blades



M2-2578




M2-2612

Scale 2:1

Art. No.	Interface	Description	Length	Shaft End	Pieces/Pkg
M2-2578	HD6	self-holding	11.5 mm	W&H	1
M2-2612	HD4	self-holding	11.5 mm	W&H	1

Twist Drills \varnothing 1.0 mm (Core Hole 1.2 Screws)

M2-3522



M2-3532

Art. No.	STERILE	Stop	Length	Shaft End	Pieces/Pkg
M2-3522	M2-3522S	5 mm	14 mm	W&H	1
M2-3532	M2-3532S	7 mm	16 mm	W&H	1

Twist Drills \varnothing 1.2 mm (Core Hole 1.5 Screws)

M2-3542




M2-3552

Art. No.	STERILE	Stop	Length	Shaft End	Pieces/Pkg
M2-3542	M2-3542S	5 mm	14 mm	W&H	1
M2-3552	M2-3552S	7 mm	16 mm	W&H	1

Twist Drills \varnothing 1.5 mm (Core Hole 1.8 and 2.0 Screws)

M2-3529



M2-3539



M2-3549



M2-3559

Art. No.	STERILE	Stop	Length	Shaft End	Pieces/Pkg
M2-3529	M2-3529S	5 mm	14 mm	W&H	1
M2-3539	M2-3539S	7 mm	16 mm	W&H	1
M2-3549	M2-3549S	9 mm	18 mm	W&H	1
M2-3559	M2-3559S	13 mm	22 mm	W&H	1

10 | 90° Screwdriver

Twist Drills Ø 1.9 mm (Core Hole for 2.3 Screws)



M2-3569



M2-3579

Art. No.	STERILE	Stop	Length	Shaft End	Pieces/Pkg
M2-3569	M2-3569S	7 mm	16 mm	W&H	1
M2-3579	M2-3579S	13 mm	22 mm	W&H	1

Twist Drills Ø 2.0 mm (Core Hole for 2.5 Screws)



M2-3589



M2-3599

Art. No.	STERILE	Stop	Length	Shaft End	Pieces/Pkg
M2-3589	M2-3589S	7 mm	16 mm	W&H	1
M2-3599	M2-3599S	13 mm	22 mm	W&H	1

Notes

DAM Medical

MANDIBLE2-09100001_v0 / © 2020-07, Medartis AG, Switzerland. All technical data subject to alteration.

MANUFACTURER & HEADQUARTERS

Medartis AG | Hochbergerstrasse 60E | 4057 Basel/Switzerland

P +41 61 633 34 34 | F +41 61 633 34 00 | www.medartis.com

SUBSIDIARIES

Australia | Austria | Brazil | France | Germany | Japan | Mexico | New Zealand | Poland | UK | USA

For detailed information regarding our subsidiaries and distributors, please visit www.medartis.com



Disclaimer: This information is intended to demonstrate the Medartis portfolio of medical devices. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. For further questions, please contact your Medartis representative (www.medartis.com). This information contains CE-marked products.
For US only: Federal law restricts this device to sale by or on the order of a physician.