

Instrumentenset voor implantologie

*Instrument kits for
Dental-Implantology*



Shown at 75% size

Surgical pliers DeBakey

TP5090

15cm, straight, TC



*DeBakey Atraumatic forceps for the safe keeping of mucosal tissue.
The fine micro-toothing prevents tissue perforation.*

TC Needle holder Castroviejo

NH5024R

18cm, straight, TC



S5081

Goldmann-Fox, 13cm
serrated, cvd.



Vestibulum Retractor Upper & Lower, 160mm

CRLA



Shown at 75% size

Periosteal Prichard

PPR3



Periosteal Molt

P9



Periodontal chisel

S13K/TG6



Bone curette Lucas

CL86



Rongeur Mini-Friedman 12cm

RF12



Double-Action Rongeur Boehler 15cm

RBOEC

Double action hinged design exerts more cutting action with each hand compression for quicker, easier bone contouring.



Instrument kit for Dental-Implantology

Shown at 75% size

Scalpel handles

10-130-5EM

Finely balanced, pen-like scalpel handle that easily rotates and maneuvers in difficult to reach areas. Scalpel blades can be used as well rights as left too.



Mirror handle with Dental mirror

M5C

Front Surface mirror # 5 - Ø 24mm.



MHE7

Dental mirror handle colour coded



Colour of silicon ring may be free choosen (grey, orange, turquoise, yellow, red, green, blue, black).

Explorer

EXS23



Periodontal probe

PCPUNC15




Colour marked in 1-mm steps 1 to 15mm

Sterilization cassette DIN size for 15 instr. + extra space

IMDIN15

290 x 184 x 34 mm

Colour of silicon rails may be free choosen (red, orange, blue, green).



Basic-Instruments kit for Sinus Lift Implantology

Shown at 75% size

Sinus Lift Instruments



4.0 mm



IMP6577SC

acute



7.4 mm



4.0 mm



IMP6578

obtuse



7.4 mm



PESLM16



IMPSC16



IMPSC26



Sterilization cassette in 1/4 DIN size for 5 instruments

IM12DIN5

72 x 184 x 34 mm

Colour of silicon rails may be free chosen
(red, orange, blue, green).



Bestelformulier ORDER FORM

Bestel-email | Order Email: info@dammedical.nl



Instrumentenset voor implantologie | Instruments kit for Dental-Implantology

Aantal Order QTY	Artikelnummer Part Code	Artikelomschrijving Article description	
	TP5090	Chirurgische Pinzette DeBakey Surgical pliers DeBakey	
	NH5024R	TC Nadelhalter Castroviejo, 18cm TC Needle holder Castroviejo, 18cm	
	S5081	Chirurgische Schere Goldman-Fox Surgical scissors Goldman-Fox	
	CRLA	Vestibulum-Abhalter für OK und UK, 160mm Vestibulum Retractor Upper & Lower, 160mm	
	PPR3	Raspatorium Prichard Periosteal Prichard	
	P9	Raspatorium Molt Periosteal Molt	
	S13K/Tg 6	Meißel für Parodontologie Periodontal chisel	
	CL86	Knochenkürette Lucas Bone curette Lucas	
	RF12	Hohlmeißelzange Mini-Friedman 12cm Rongeur Mini-Friedman 12cm	
	RBOeC	Double-Action Hohlmeißelzange Boehler 15cm Double-Action Rongeur Boehler 15cm	
	10-130-5em	Skalpelgriff, rund, gerade Scalpell handle, round, straight	
	mHe7 + m5C	Spiegelgriff + Front Surface Dentalspiegel # 5 Mirror handle + Front Surface dental mirror # 5	
	eXS23	Zahnsonde Explorer	
	PCPUNC15	Parodontometer UNC 15 Periodontal probe UNC 15	
	ImDIN15	Sterilisationskassette, DIN, 15 Instr. + Extraplatz Sterilization cassette DIN, 15 instr. + extra space	
	KIT-IMP1	Complete kit Complete Kit as above listed	

Basis-Instrumentenkit voor sinus lift | Basic-Instruments kit for Sinus Lift Implantology

Aantal Order QTY	Artikelnummer Part Code	Artikelomschrijving Article description	
	ImP6577SC	Sinus Lift instrument	
	ImP6578	Sinus Lift instrument	
	PeSIm16	Sinus Lift instrument	
	ImPSC16	Sinus Lift instrument	
	ImPSC26	Sinus Lift instrument	
	Im14DIN5	Sterilisationskassette ¼ DIN für 5 Instrumente Sterilization cassette ¼ DIN for 5 instruments	
	KIT-IMP5	Complete kit Complete Kit as above listed	

Fachhändler-/ Depot-Anschrift eintragen (wenn vorhanden)
Dealer company address (if available)

Praxis | Clinic name

Ansprechpartner | Contact person

Straße | Address

PLZ | ZIP

Ort | City

Telefon | Phone

Datum | Order date

Unterschrift | Signature