



Less is more

Discover our minimally invasive solutions with the family of extra-short implants, starting at just 4.5 mm in total length.

All the **SCIENCE**,
all the **EXPERIENCE**,
all the **RELIABILITY**,
all the **BENEFITS**,



now in
4,5mm



extra short implants

Less is more




Atraumatic solution for severe atrophies in both jaws.

Simple and fast implant site preparation.

Predictable surgery, with no need for complex and invasive techniques.

The key is in the ***Biomechanics***

Biomechanical studies back the use of short implants, in comparison with longer implants, proving the capacity to dissipate the generated tensions both from vertical and lateral loads, in the first millimeters of the neck.



4.5 mm length implants are available on the
CORE Narrow platform and the Universal Plus platform

Narrow/CORE Platform

Ø3,5 mm platform diameter

Diameters: 3,3 | 3,5 | 3,75 | 4
4,25 | 4,75 mm




Universal Plus Platform

Ø4,1 mm platform diameter

Diameters: 4,5 | 5 | 5,5 | 6 mm

**What prosthetic solution
is the right one for
these implants?**



**Only multiple restorations screw-retained to
the specific Multi-Im transepithelials*, both
straight and expanded, to obtain a different
prosthetic emergence.**



*** Important: The use of a non-specific component will mean a bad fit between pieces.
The gap created between components will break the hermetic seal and passive fit.**

Clinical case

*¿Do you want to see
a clinical case of the
4.5 mm length implant?*

We invite you to see it



BTI Biotechnology Institute

info@bti-biotechnologyinstitute.com

bti-biotechnologyinstitute.com

